

UNIVAC DR. FRENCH PORCELAIN POSTERIORIS

OCCLUSAL FORM


- Poly-Plane
- Designed to reduce lateral thrust created by masticating forces

MOLDS

- 39 - long medium, and short
- 42 - long, medium, and short
- 45 - medium and short
- 48 - medium

COLORS

Univac Porcelain B1, Y1, Y2, Y3, Y4, Y5, Y6



This mold series for all posteriors include:

- A. Mold 39 = 28 millimeters
- B. Mold 42 = 30 millimeters
- C. Mold 45 = 32 millimeters
- D. Mold 48 = 34 millimeters

Note: The approximate width of the four posteriors on one side (2 molars and 2 bicuspids)

UNIVAC DR. FRENCH PORCELAIN POSTERIORIS

Designed and developed over 8 years, Dr. French's Posteriors are designed with the following features:

They resemble the natural in so far as function will permit. * Uppers are made so they accept the movements of the lowers and function with them. * Balanced occlusion can be easily obtained without interference. If we have contact on both sides of the mouth during working excursion, we distribute the force of the bite more evenly over the entire area instead of concentrating it on one point. * They are easy running and won't lock or trip. * They allow for normal settling without disturbing articulation. * They are efficient in mastication. * They meet mechanical requirements. That is, they function under edentulous conditions, in accordance with sound mechanical laws.

