



1-800-661-2044

Optimum dental esthetics in porcelain teeth

UNIVAC® Teeth

- ▼ Universal Dental mold system simplifies articulation: mold numbers correspond to millimeter width of 6 upper anteriors carded flat
- ▼ Dentin-enamel color blending for true-to-life esthetics
- ▼ Also available in POLYCHROME® characterized anterior teeth

ANTERIORS		POSTERIORS			WIDTH OF 6 ANTERIORS CARDED FLAT		UPPER CENTRALS	
UPPER	LOWER	BIO-MECH.	NIC.	DR. FRENCH'S	UPPER	LOWER	WIDTH	LENGTH EXC. OF COLLAR

C - DOMINANTLY CURVED LABIALS

L39/C	L39	M39	L39	L39	39.2	30.0	7.2	9.6
M39/C	M39	M39	M39	M39	39.2	30.0	7.2	8.5
S39/C	S39	S39	S39	S39	39.2	30.0	7.2	7.6
L42/C	L42	M42	L42	L42	42.2	32.0	7.7	10.3
M42/C	M42	M42	M42	M42	42.2	32.0	7.7	9.1
S42/C	S42	S42	S42	S42	42.2	32.0	7.7	8.1
L45/C	L45	M45	L45	L45	45.2	34.0	8.2	10.9
M45/C	M45	M45	M45	M45	45.2	34.0	8.2	9.6
S45/C	S45	S45	S45	S45	45.2	34.0	8.2	8.6
L48/C	M48	M48	L48	L48	48.2	36.0	8.7	11.6
M48/C	M48	M48	M48	M48	48.2	36.0	8.7	10.2
S48/C	S48	S48	S48	S48	48.2	36.0	8.7	9.2

FA - DOMINANTLY FLAT LABIALS

L39/FA	L39	M39	L39	L39	39.6		7.2	9.8
M39/FA	M39	M39	M39	M39	39.6		7.2	9.0
S39/FA	S39	S39	S39	S39	39.6		7.2	8.0
L42/FA	L42	M42	L42	L42	42.6		7.7	10.7
M42/FA	M42	M42	M42	M42	42.6		7.7	9.4
S42/FA	S42	S42	S42	S42	42.6		7.7	8.8
L45/FA	L45	M45	L45	L45	45.6		8.2	11.1
M45/FA	M45	M45	M45	M45	45.6		8.2	10.0
S45/FA	S45	S45	S45	S45	45.6		8.2	9.6
L48/FA	M48	M48	L48	L48	48.6		8.8	11.6
M48/FA	M48	M48	M48	M48	48.6		8.8	10.3
S48/FA	S48	S48	S48	S48	48.6		8.8	9.6

SUPPLEMENTAL MOLDS

M42/CS	M39	M42	M42	M42	42.4		8.0	10.0
L44/FS	M42	M45	M45	M45	43.8		8.2	10.7
M44/FS	M42	M42	M42	M42	44.4		8.0	9.7
L46/FS	M42	M45	M45	M45	45.0		8.5	11.2
L46/CS	M42	M45	M45	M45	45.4		8.5	10.8

UNIVAC® AND UNIVAC POLYCHROME® PORCELAIN ANTERIORS

C DOMINANTLY CURVED LABIALS

All measurements are in millimeters and approximate. The lengths of the central incisors are exclusive of the collar.

LONG

MEDIUM

SHORT

39 mm			Width 7.2, Length 9.5			Width 7.2, Length 8.5			Width 7.2, Length 7.6
42 mm			Width 7.7, Length 10.3			Width 7.7, Length 9.1			Width 7.7, Length 8.1
45 mm			Width 8.2, Length 10.9			Width 8.2, Length 9.6			Width 8.2, Length 8.6
48 mm			Width 8.7, Length 11.6			Width 8.7, Length 10.2			Width 8.7, Length 9.2

FA DOMINANTLY FLAT LABIALS

All measurements are in millimeters and approximate. The lengths of the central incisors are exclusive of the collar.

39 mm			Width 7.2 Length 9.8			Width 7.2, Length 9.0			Width 7.2, Length 8.0
42 mm			Width 7.7, Length 10.7			Width 7.7, Length 9.4			Width 7.7, Length 8.8
45 mm			Width 8.2, Length 11.1			Width 8.2, Length 10.0			Width 8.2, Length 9.6
48 mm			Width 8.8, Length 11.6			Width 8.8, Length 10.3			Width 8.8, Length 9.6

The Coordinate Size Mold System

All of the anterior molds are arranged systematically in an orderly plan that classifies each set by:

- Dominant labial character C and FA
- Length - long, medium or short
- ▲ The width of six (6) anteriors, carded flat, is in millimeters.

SUPPLEMENTAL MOLDS

All measurements are in millimeters and approximate. The lengths of the central incisors are exclusive of the collar.



AVAILABLE IN LONG AND MEDIUM VERTICALS WITH FLAT (FS) AND CURVED (CS) LABIAL SURFACES

LOWERS

All measurements are in millimeters and approximate. The lengths of the central incisors are exclusive of the collar.




UNIVAC PORCELAIN POSTERIOR

UNIVAC N.I.C. (non-interfering cusps) PORCELAIN POSTERIOR

M39



M45



M42




M48



UNIVAC BIO-MECHANICAL PORCELAIN POSTERIOR

M39




M45



M42



M48



UNIVAC N.I.C. (non-interfering cusps) PORCELAIN POSTERIOR

OCCLUSAL FORM

- Semi-anatomical, 20° posterior
- Designed to avoid anteroposterior lock, with a narrow bucco-lingual dimension of the lower molds

COLORS

Univac Porcelain B1, Y1, Y2, Y3, Y4.

MOLDS

39 - long, medium, and short
42 - long, medium, and short
45 - long, medium, and short
48 - long, medium, and short

The mold series for all posteriors include:

- A. Mold 39 = 28 millimeters
- B. Mold 42 = 30 millimeters
- C. Mold 45 = 32 millimeters
- D. Mold 48 = 34 millimeters

UNIVAC BIO-MECHANICAL PORCELAIN POSTERIOR

OCCLUSAL FORM

- Flat Plane, 0° posteriors
- Biologic in function, mechanical in design. May be set on a curved or flat occlusal plane.

COLORS

Univac Porcelain B1, Y1, Y2, Y3, Y4.

MOLDS

39 - Medium and short
42 - Medium and short
45 - Medium and short
48 - Medium and short



Note: Short molds are designed for a minimum of grinding.

Note: The approximate width of the four posteriors on one side (2 molars and 2 bicuspids)